

# Product Specifications

## Tissue-Tek® Base Mold System



Base Mold System

### Product description

Item Code	Product Description	Sales quantity
4067	Torpedo Base Mold for 6-Chamber Cassette (for item 4073 - 4077)	12/pcs
4121	Base Mold 16 x 16 x 6 mm	12/pcs
4122	Base Mold 22 x 22 x 6 mm	12/pcs
4123	Base Mold 32 x 25 x 6 mm	12/pcs
4124	Base Mold 38 x 25 x 6 mm	12/pcs
4131	Base Mold 22 x 22 x 12 mm	12/pcs
4132	Base Mold 32 x 25 x 12 mm	12/pcs
4133	Base Mold 38 x 25 x 12 mm	12/pcs
4161	Base Mold 7 x 7 x 5 mm	12/pcs
4162	Base Mold 15 x 15 x 5 mm	12/pcs
4163	Base Mold 24 x 24 x 5 mm	12/pcs
4164	Base Mold 30 x 24 x 5 mm	12/pcs
4165	Base Mold 37 x 24 x 5 mm	12/pcs
4166	Mega-Cassette® Base Mold (for item 4173)	6/pcs
7840	Super Mega-Cassette® Base Mold 36x36x10 mm,	5/pcs
7841	Super Mega-Cassette® Base Mold 60x45x15 mm,	5/pcs
7055	Paraform® Standard Base Mold	12/pcs
7056	Paraform® Biopsy/Orientation Base Mold	12/pcs
7057	Paraform® Biopsy 13 x 13 Base Mold	12/pcs

Tissue Embedding

Tissue-Tek® Base Mold System



**Product Specifications**  
**Tissue-Tek® Super Mega Cassette**

**Product specifications**

<b>1. Application</b>	
	Designed to embed specimen into a paraffin block. Different sizes of basemolds can be used for different tissue dimensions.

<b>2. Related Items</b>	
Item Code	Product Description
VSP-5-9060-9067	Super Mega-Cassette® System
4173	Mega-Cassette®
4073-4077	Biopsy 6-Chamber Cassette
4306-4314	Mesh Biopsy Cassette
7019-7024	Paraform® Cassette System
	Tissue-Tek® III (Biopsy) Uni-Cassette® System
	Tissue-Tek® II Processing/Embedding Cassette System

<b>3. Regulatory information</b>	
Intended Use	Designed to embed specimen into a paraffin block.
CE status	CE
MSDS available	<b>Yes</b> <input type="checkbox"/> (xx Languages) <b>No</b> <input checked="" type="checkbox"/>

<b>4. Packaging</b>	
Package description	Carton box
Sales quantity	See product description
Dimensions (LxWxH)	<div style="display: flex; flex-direction: column; align-items: flex-start;"> <div style="margin-bottom: 20px;"> <p><b>4067</b></p> <p>X: 35 mm Y: 7 mm Z: 5 mm</p>  </div> <div> <p><b>4121</b></p> <p>X: 16mm Y: 16mm Z: 6mm</p>  </div> </div>

Product Specifications  
Tissue-Tek® Super Mega Cassette

**4122**

X: 22mm  
Y: 22mm  
Z: 6mm



**4123**

X: 32mm  
Y: 25mm  
Z: 6mm



**4131**

X: 22mm  
Y: 22mm  
Z: 12mm



Product Specifications  
Tissue-Tek® Super Mega Cassette

**4132**

X: 32mm  
Y: 25mm  
Z: 12mm



**4133**

X: 38mm  
Y: 25mm  
Z: 12mm



**4161**

X: 7mm  
Y: 7mm  
Z: 5mm



Product Specifications  
Tissue-Tek® Super Mega Cassette

**4162**

X: 15mm  
Y: 15mm  
Z: 5mm



**4163**

X: 24mm  
Y: 24mm  
Z: 5mm



**4164**

X: 30mm  
Y: 24mm  
Z: 5mm



11

Product Specifications  
Tissue-Tek® Super Mega Cassette

**4165**

X: 37mm  
Y: 24mm  
Z: 5mm



12

**4166**

X: 30 mm  
Y: 22 mm  
Z: 12 mm



**7055**

X: 32mm  
Y: 24mm  
Z: 5mm



Product Specifications  
Tissue-Tek® Super Mega Cassette

**7056**

X: 25mm  
Y: 18mm  
Z: 5mm



**7057**

X: 13mm  
Y: 13mm  
Z: 5mm



**7840**

X: 36mm  
Y: 36mm  
Z: 10mm




**Product Specifications**  
**Tissue-Tek® Super Mega Cassette**

<b>7841</b>	<p>X: 65mm  Y: 45mm  Z: 15mm</p>
Package dimensions (LxWxH)	100 x 60 x 35 mm
Shipping dimensions (LxWxH)	100 x 60 x 35 mm
Shipping weight (approx.)	0.15 kg

**5. Labeling**

HIBC code	+H905 <b>Item Code</b> 47
Statistical code	73102990



*Label item 7055*

**6. Storage and transport**

Storage conditions:	
• Temperature	NA
• Humidity	NA
• Sunlight	NA
Transport conditions:	
• Temperature	NA
• Humidity	NA
• Sunlight	NA
Handle with care (fragile)	<b>Yes</b> <input type="checkbox"/> <b>No</b> <input checked="" type="checkbox"/>
Dangerous Good	<b>Yes</b> <input type="checkbox"/> <b>No</b> <input checked="" type="checkbox"/>
Shelf life	5 years

**7. Material specifications**

Components	Stainless steel
Chemical resistance	Reacting with acid
Heat resistance	NA
Cleaning	Xylene/ Tissue-Clear® or other xylene substitutes Alcohol, isopropyl alcohol.



**Product Specifications**  
**Tissue-Tek® Super Mega Cassette**

Latex free	<b>Yes*</b> <input checked="" type="checkbox"/> <b>No</b> <input type="checkbox"/> *only when official statement is available
------------	--

<b>8. Terms &amp; Conditions</b>	
Warranty	6 months after receiving by buyer
Terms of Delivery	F.O.B. Warehouse, The Netherlands

<b>9. Additional information</b>	