

Product Specifications

Tissue-Tek® FormaGO



Tissue-Tek® FormaGO Formalin Absorption Granules

Product description

Item Code	Product Description	Sales quantity
9158	Tissue-Tek® FormaGO Formalin Absorption Granules	1x2 (2x900gr.)

Product specifications

1. Application	
	The granules are ideal for large formalin spills, and will protect laboratory personnel from exposure to dangerous vapors. The light weight provides convenient handling. After distributing the granules over the spill, they will absorb the spilled formalin. The granules retain their form, making it easier to pick them up with a broom or brush. Use gloves as the granules can cause skin irritation.

2. Related Instrument	
Item Code	Product Description

3. Regulatory information	
Intended Use	The Formalin Absorption Granules are used to clean up larger formalin spills in the laboratory. Granules that are put on the spill will absorb the formalin, making it easy to clean up.
CE status	
MSDS available	Yes <input checked="" type="checkbox"/> No <input type="checkbox"/>


Tissue Processing

Tissue-Tek® FormaGO

Product Specifications
Tissue-Tek® FormaGo

4. Packaging	
Package description	Pot of 900gr.
Sales quantity	1x2
Dimensions (LxWxH) or (ØxH)	10x10x21cm
Package dimensions (LxWxH)	
Shipping dimensions (LxWxH)	
Shipping weight	2,27 kg

5. Labeling	
HIBC code	
EAN code	
Statistical code	



Label 9158

6. Storage and transport	
Storage conditions: <ul style="list-style-type: none"> • Temperature • Humidity • Sunlight 	Store in a closed container at room temperature, 15-30°C.
Transport conditions: <ul style="list-style-type: none"> • Temperature • Humidity • Sunlight 	Not regulated. Excessive cold/ heat/ light should be avoided.
Handle with care (fragile)	Yes <input checked="" type="checkbox"/> No <input checked="" type="checkbox"/>
Dangerous Good	Yes <input type="checkbox"/> No <input type="checkbox"/>
Shelf life	

7. Material specifications	
Components	Mixture of trisodium phosphate, sodium metabisulfite and diatomaceous earth.
Chemical resistance	Contact with acids will release sulfur dioxide
Heat resistance	Excessive cold/ heat should be avoided
Cleaning	Clean area with water after use
Latex free	Yes* <input checked="" type="checkbox"/> No <input type="checkbox"/> *only when official statement is available

Product Specifications
Tissue-Tek® FormaGo

8. Terms & Conditions	
Warranty	6 months after receiving by buyer
Terms of Delivery	F.O.B. Warehouse, The Netherlands

9. Additional information	