

---

**Tissue-Tek®**

Microscope Slides portfolio

*Imagine... a suitable  
slide for every situation*



continuous innovation for pathology



# A SUITABLE MICROSCOPE SLIDE FOR EVERY SITUATION.

For everyday use in routine Pathology, for use in IHC/ ISH, in case of special stains or other situations where tissue is known to fall of slides easier; Sakura has the solution.

## FULL MICROSCOPE SLIDES PORTFOLIO

Tissue-Tek® Microscope Slides are available as:

- Non-Adhesive
- Adhesive
- Advanced Adhesive

Different slides are available for:

- Label printers
- Inkjet printers
- Thermal printers

### Available finishing of slide edges:

- **Cut edges** – economical solution for routine applications
- **Ground edges 90°** – the ground edge finish serves as a safety precaution to minimize the chance of cuts
- **Ground edges 90° in combination with clipped corners** – is most suitable for Sakura automatic equipment; The smooth ground edges make slides move easier through instruments, while the clipped corners give less chance of slides getting stuck while moving



## SMARTWRITE™

SMARTWrite™ Non-Adhesive Microscope Slides

SMARTWrite™ Adhesive Microscope Slides

Slides specifically designed to be used in combination with thermal printers and inkjet printers, or use pencil writing.

- Available with white colored tab only
- Available with ground edges and clipped corners

## GLAS™

Non-Adhesive Microscope Slides

- For everyday use in routine Pathology
- Available with cut edges or ground edges
- Available with or without frosted ends
- Suitable for use with label printers
- Straightforward slides without colored tabs or spacers



## AUTOWRITE®

**AutoWrite® Non-Adhesive Microscope Slides**

**AutoWrite® Adhesive Microscope Slides**

**AutoWrite® Advanced Adhesive Microscope Slides**

Slides specifically designed to be used in combination with Inkjet printers.

- Silane coated
- Available with 7 different colored tabs
- Available with ground edges and clipped corners
- Extremely suitable for IHC and ISH

## DEDICATED TO HISTOPATHOLOGY

Sakura Finetek improves the laboratories. By offering a total solution for cryotomy, Sakura Finetek makes your laboratories ready to deliver standardized, superior quality, results. As the innovative company in histopathology, Sakura Finetek is continuously looking for possibilities to improve laboratories... and succeeds in offering the only complete package to meet present and future standards in the histopathology laboratories around the world.

## PRODUCT RANGE

### Tissue-Tek® AutoWrite® Non-Adhesive

9562	Ground Edges 90° Clipped Corner White
9563	Ground Edges 90° Clipped Corner Blue
9564	Ground Edges 90° Clipped Corner Pink
9565	Ground Edges 90° Clipped Corner Yellow
9566	Ground Edges 90° Clipped Corner Green
9567	Ground Edges 90° Clipped Corner Orange
9568	Ground Edges 90° Clipped Corner Purple

### Tissue-Tek® AutoWrite® Advanced Adhesive

9618	Ground Edges 90° Clipped Corner White
------	---------------------------------------

### Tissue-Tek® AutoWrite® Adhesive

9597	Ground Edges 90° Clipped Corner White
9598	Ground Edges 90° Clipped Corner Blue
9599	Ground Edges 90° Clipped Corner Pink
9600	Ground Edges 90° Clipped Corner Yellow
9601	Ground Edges 90° Clipped Corner Green
9602	Ground Edges 90° Clipped Corner Orange
9603	Ground Edges 90° Clipped Corner Purple



## PRODUCT RANGE

### Tissue-Tek® Glas™

9531	Glas™ Non-Adhesive Cut Edge
9532	Glas™ Non-Adhesive Ground Edges 90
9533	Glas™ Non-Adhesive Ground Edges 90° Frosted
9534	Glas™ Non-Adhesive Cut Edges Frosted
9579	Glas™ Coverslips 24 x 55 mm
9580	Glas™ Coverslips 24 x 40 mm
9581	Glas™ Coverslips 24 x 50 mm
9582	Glas™ Coverslips 24 x 60 mm

### Tissue-Tek® Glas™

9641	Glas™ Super Mega Slide Cut Edges 76 x 52
9642	Glas™ Super Mega Coverslips 60 x 50 mm
9643	Glas™ Super Mega Coverslips 76 x 52 mm

### Tissue-Tek® SMARTWrite™

9585	Non-Adhesive Ground Edges 90° Clipped Corner White
9721	Adhesive Ground Edges 90° Clipped Corner White

**Contact your Sakura Sales Representative to order your Sakura microscope slides.**

Please visit our website [sakura.eu](http://sakura.eu)

Sakura Finetek Europe B.V., Flemingweg 10A, 2408 AV Alphen aan den Rijn,  
P.O. Box 362, 2400 AJ Alphen aan den Rijn, The Netherlands

