

---

## Tissue-Tek® VIP® 6 AI

Vacuum Infiltration Processor

*Imagine... outstanding  
reliability meeting  
exceptional safety and  
ultimate quality*



continuous innovation for pathology



## PRODUCT SPECIFICATIONS

<b>Capacity</b>	max. 300 cassettes
<b>Dimensions (W x D x H)</b>	61 x 68 x 133 cm
<b>Weight</b>	175 kg
<b>Reagent reservoirs</b>	13 (including waste and cleaning)
<b>Paraffin reservoirs</b>	4
<b>Waste containers</b>	2 (1 reagent and 1 paraffin wax)
<b>Bulk containers</b>	2 (10,8 L each)
<b>Reagent Management System</b>	3 reagents including paraffin wax
<b>Level sensors</b>	4 (150-cassettes, 300-cassettes, cleaning and overflow)
<b>Fume control</b>	Fume condenser system and 2 disposable activated carbon filters

<b>Screen</b>	Colour touch screen
<b>Program memory</b>	50 processing programs, 3 user-defined cleaning programs
<b>Password protection</b>	2 levels, 4 managers, and 20 users



## ITEM DESCRIPTIONS

<b>6040</b>	VIP® 6 AI Vacuum Infiltration Processor; 100V, 50/60Hz	<b>6036</b>	1 set of bottle labels, 10 each
<b>6041</b>	VIP® 6 AI Vacuum Infiltration Processor; 115V, 50/60Hz	<b>6045</b>	1 Set of paraffin waste bags
<b>6042</b>	VIP® 6 AI Vacuum Infiltration Processor; 230V, 50/60Hz	<b>6160</b>	Activated carbon filter
<b>1550</b>	Scraper, large	<b>7109</b>	Basket transportation tray
<b>6034</b>	Scraper, small	<b>Optional items</b>	
<b>5267</b>	150- cassette baskets with dividers	<b>5265</b>	150- cassette basket with dividers and handles
<b>5269</b>	Basket handles for 150- cassette basket (5267)	<b>5266</b>	75- cassette basket with dividers and handles
<b>6035</b>	Reagent bottles, complete, 1/ EA	<b>5268</b>	75- cassette baskets with dividers
		<b>5270</b>	Basket handles for 75- cassette basket (5268)